

# Diagnostics for SXFEL

---

- diagnostics specific for single pass FEL
- standard diagnostics tools for e-bunch measurement
- coherent radiation measurement of short bunches
- single-shot e-bunch measurement
- outlook---Prof. Yen-Chieh Huang

**Bian Yu**  
2015/9/9

- no photon diagnostics
- personal perspective



Most of the materials are available on the Internet, & this ppt is only for internal communication in FEL division at SINAP.

# Diagnostics specific for single pass FEL

Main parameters of SXFEL & SDUV		
Parameters	SXFEL	SDUV
Electron Energy	0.84 GeV	130MeV
Normalized emittance	< 2.5mm.mrad	< 6.0mm.mrad
Slice energy spread	< 0.02%	\
Peak current/ bunch length/charge	500 A /1ps/ 500pC	100A/100fs~3ps
Wavelength of laser seed	~265nm	1047nm/800nm
Scheme	HGHG	HGHG/EEHG
Cascading scheme	265-44-8.8nm;	\
FEL peak power	>100MW	~MW
FEL wavelength	8.8 nm	200~1200nm
FEL pulse length	100~150 fs (FWHM)	8ps/~130fs/~100fs

FEL pulse length and e-bunch length

# Diagnostics specific for single pass FEL




FEL power depends exponentially on beam parameters (peak current, emittance ... )

Measure, control and stabilize beam parameters such that optimum FEL performance is achieved


# Diagnostics specific for single pass FEL

Bunch compression for high peak currents has non-linear components  $\Rightarrow$  complex phase space distributions




$\Delta E$  [MeV]

Expected long. bunch shape at LCLS, 'double horn' due to wake fields







$I$  [A]



Only fraction of the total charge will 'lase', diagnostic has to be sensitive to this fraction

# Diagnostics specific for single pass FEL



ideal diagnostics

- ultimate resolution
- comprehensive
- immediate feedback
- non-destructive

Projected parameters are of limited use !


Diagnostics has to reveal details of the bunch structure

slice emittance, bunch profile, slice energy spread, bunch position

# Standard diagnostics tools for e-bunch measurement

## Transverse deflecting cavities (TCAV)

- Adds z-position dependent transverse kick to bunch
- Phase advance to screen vertical streak of longitudinal bunch structure



- single bunch capable
- destructive (sacrifice 1 bunch)
- slow read out (imaging)


Typical Resolution: 20-50fs

Resolution depends on cavity power, beam energy and machine optics

# Standard diagnostics tools for e-bunch measurement

## Electro -Optic (EO) Techniques

Intra-beamline measurement of the bunch Coulomb field



- Field induced refractive index change
- Polarization-modulation of probing laser
- Temporal structure of Coulomb field impressed to ellipticity of optical pulse

- scanning techniques
- **single-shot techniques**

Limitations:

- high frequency cut-off due to finite distance to beam
- velocity mismatch of FIR and optical propagation in EO crystal
- phonon resonances of EO material

# Standard diagnostics tools for e-bunch measurement

Make the electrons radiate coherently ...

*spectral energy density*

source characteristics (CSR, CTR, CDR...)

$$\frac{dU}{d\omega} = CN^2 \left| F_{long}(\omega) \right|^2 T(\omega, \gamma, r_b, \theta, source)$$

$$F_{long}(\omega) = \int_{-\infty}^{\infty} \tilde{\rho}(t) \exp(-i\omega t) dt$$

normalized charge density

-integral intensity



compression factor, effective bunch length

-spectral resolved intensity



bunch structure, longitudinal beam profile



## CTR measurement of short bunches----in frequency domain



### Single-particle transition radiation (TR)

$$I_0(\omega) = \frac{d^2 U_{1e}^-}{d\Omega dk} = \frac{r_e mc}{\pi^2} \frac{\beta^2 \sin^2 \theta}{(1 - \beta^2 \cos^2 \theta)^2} \approx \frac{e^2}{\pi^2 c} \frac{\theta^2}{(\gamma^{-2} + \theta^2)^2}$$

Frequency independent




Ideal for analysis




For SDUV, e-bunch energy is 130MeV.

# Coherent transition radiation(CTR) from e-bunch



Bunch length,  $\sigma_z$

Incoherent radiation  $\lambda < \sigma, I \propto N$



Coherent radiation  $\lambda > \sigma, I \propto N^2$

Far-field approximation:

$$I_{coherent}(\omega, \sigma_z) = I_0(\omega) [N + \underline{N(N-1)} F_b(\omega, \sigma_z)] \\ \sim N^2 I_0(\omega) F_b(\omega, \sigma_z)$$

$I_{coherent}(\omega, \sigma_z)$  : Power of coherent radiation

$N$ : Number of particle in a bunch

$F_b(\omega, \sigma_z)$  : Bunch form factor (BFF)

$I_0(\omega)$ : Radiation emitted from an electron

$\sigma_z$  : Bunch length

Bunch form factor :

$$F_b(\omega) = |f(\omega)|^2$$



$$f(\omega) = \sum_{j=1}^N e^{i\omega t_j} \cong N \int_{-\infty}^{\infty} \rho(t) e^{i\omega t} dt$$

**stack in same phase!**

# Coherent transition radiation(CTR) from e-bunch

For a gaussian beam the longitudinal beam profile is:

$$\rho(t) = \frac{1}{\sqrt{2\pi}\sigma} \exp(-t^2 / 2\sigma^2) \quad \xrightarrow{\text{BFF}} \quad F_b(\omega) = \exp(-\sigma^2 \omega^2)$$




short e- bunch  $\Leftrightarrow$  radiation at higher frequency


# Facility in Osaka University (<20fs)

## Bunch length measurement system (in a vacuum.)


M: mirror, OAP: off-axis parabolic mirror,  
BS: beam splitter (high resistivity silicon or KBr), MCT: HgCdTe detector



Michelson interferometer



Infrared detectors




Pictures are from Professor Yang Jinfeng in Osaka University

Two different infrared detectors are used!  
(Bolometer: low-freq. radiation, MCT: higher freq. radiation)


# e-bunch length estimation

Exp. Interferogram data,  $i_{eff}(t)$




FT

Freq. spectrum,  $I_{eff}(\omega)$





$S(\omega)$



assuming Gaussian distr.

$$\sigma = 8.9 \text{ fs}$$



# Longitudinal beam profile reconstruction

Using Kramers-Kronig (KK) transform

$$F_b(\omega) = I_{eff}(\omega) / S(\omega)$$

$$\phi(\omega_1) = -\frac{\omega_1}{\pi} \int_0^\infty \frac{\ln[F_b(\omega)/F_b(\omega_1)]}{\omega^2 - \omega_1^2} d\omega$$

$$\rho(t) = \int_{-\infty}^\infty \sqrt{F_b(\omega)} \exp[-i\{\omega t - \phi(\omega)\}] d\omega$$


$\rho(t)$  : Charge distribution

$\phi(\omega)$  : Phase at each frequency

$I_{eff}(\omega)$  : Effective frequency spectrum

$S(\omega)$  : Sensitivity response function

$F_b(\omega)$  : Bunch form factor



# Longitudinal beam profile reconstruction

Using Kramers-Kronig (KK) transform

$$F_b(\omega) = I_{eff}(\omega) / S(\omega)$$

$$\phi(\omega_1) = -\frac{\omega_1}{\pi} \int_0^\infty \frac{\ln[F_b(\omega)/F_b(\omega_1)]}{\omega^2 - \omega_1^2} d\omega$$

$$\rho(t) = \int_{-\infty}^\infty \sqrt{F_b(\omega)} \exp[-i\{\omega t - \phi(\omega)\}] d\omega$$


$\rho(t)$  : Charge distribution

$\phi(\omega)$  : Phase at each frequency

$I_{eff}(\omega)$  : Effective frequency spectrum

$S(\omega)$  : Sensitivity response function

$F_b(\omega)$  : Bunch form factor



# Longitudinal beam profile reconstruction

Using Kramers-Kronig (KK) transform

$$F_b(\omega) = I_{eff}(\omega) / S(\omega)$$

$$\phi(\omega_1) = -\frac{\omega_1}{\pi} \int_0^\infty \frac{\ln[F_b(\omega)/F_b(\omega_1)]}{\omega^2 - \omega_1^2} d\omega$$

$$\rho(t) = \int_{-\infty}^\infty \sqrt{F_b(\omega)} \exp[-i\{\omega t - \phi(\omega)\}] d\omega$$


$\rho(t)$  : Charge distribution

$\phi(\omega)$  : Phase at each frequency

$I_{eff}(\omega)$  : Effective frequency spectrum


$S(\omega)$  : Sensitivity response function

$F_b(\omega)$  : Bunch form factor





# CSR measurement of short e-bunch



Depending on compression scheme

Coherent effects create spectral substructure

Micro-bunching can produce  $\sim$  few  $\mu\text{m}$  coherent radiation

Technical implications

- diamond windows to accelerator vacuum
- no radiation transport in (humid) air
- broad wavelength range to cover, **SINGLE SHOT**



# single-shot e-bunch measurement

Classical : Michelson type interferometers

- scanning devices, no single shot
- complex unfolding procedure (autocorrelation function)


Single shot spectrometers:

dispersive elements & multichannel detector



# single-shot e-bunch measurement

## Single shot multichannel detectors



- + commercial
- + fast
- + sensitive
- cryogenic device
- very expensive

## Requirements :

- fast, 200ns for XFEL bunch spacing
- uniform spectral response
- broadband ( $1\ \mu\text{m}$  -1mm)




## Pyro-electric line detector


- + room temperature
- + sensitivity  $\sim 300\ \text{pJ}$

# single-shot e-bunch measurement


## Smith-Purcell radiation measurements



Measurement at 45 MeV, FELIX




see PRST 9,092801 (2006)



Results of a run at 28.5 GeV from SLAC are currently being analyzed.

# outlook---Prof. Yen-Chieh Huang

## 'FROG' + Pyro-electric line detector + Reference Laser gun




### Advantages

- room-temperature
- inexpensive
- non-THz wave

A sketch of **F**requency-**r**esolved **o**ptical **g**ating system



# Few femtosecond, few kiloampere electron bunch produced by a laser-plasma accelerator



Thank you for your attention !